MedLinket

Pulse Oximeter BZ1321B REV:A.1

Welcome Guide

Model: AM801 (REF: AM-601)

How to begin?

1.Open battery compartment cover and insert two new AAA alkaline batteries as indicated below to ensure the polarity (<+> and <->) of the batteries are correct and then close it. 2*AAA alkaline battery

2.Open the clip of the pulse oximeter, insert a finger (The preferred application is middle finger and index finger), as shown below.



3. Make sure the sensor LED is above the fingernail as shown in the below.



Content Foreword Warning Latex Content Statement About This Manual Section1-Overview Intended Use .. About the Pulse Oximeter Identification of Front Panel and Symbols2 Equipment Symbols Technical Specifications -------3 Product parts and accessories Principle of Measurement Clinical Restrictions ------5 Attentions ... Unpacking and Inspection Included in the package Installation and Verification ----- 6 General Operation------ 7 Safety Function Setting Introduction Menu Setting...q Section 3- Troubleshooting 11 Section 4- Electromagnetic Environment Electromagnetic Environment 12 Section 5- Measurement Validation 13 Subject Demographics 13 ARMS Results 13 Graphs .. Section 6- Service and Maintenance ----- 15 Cleaning and Disinfecting 15 Calibrating Repairing and Maintenance 15 Disposal 16 Warranty..... Qualified certificate Foreword The Pulse Oximeter manual is intended to provide information for proper operation and maintenance. General knowledge of monitoring and understanding of the features and functions of the Pulse Oximeter Monitor are prerequisites for proper use. Please read these instructions carefully before using this equipment. The manual describing the operating procedures should be followed strictly. Failure to follow these instructions can cause measuring abnormality, equipment damage and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The Pulse Oximeter is a medical device, and can be used repeatedly. Warning / Warnings are identified by the WARNING symbol shown above. • Explosion hazard. Do not use the PULSE OXIMETER in the presence of

data values but these data should not be used as a basis for a clinical diagnosis. • Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.

flammable anesthetics mixed with air, or with oxygen, or nitrous oxide.

skin integrity changes, move the sensor to another site.

Do not spray, pour, or spill any liquid on the PULSE OXIMETER, its accessories,

skin condition affects the ability of the skin to tolerate sensor placement, it may

be necessary to change the sensor site more frequently with some patients. If

prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors

operate without risk of exceeding 41 °C on the skin if the initial skin temperature

Be aware that following removal of the sensor from the patient, it is possible that

• Sensors must be moved to a new site at least every 4 hours. Because individual

At elevated ambient temperatures, patient skin could be severely burned after

Equipment Symbols

Equipmen	Equipment Symbols		
Â	Caution	86kpa	Atmospheric pressure limitation (Storage, Transportation and Operating)
NON	Non-sterile	×	Type BF applied part
Ĩ	Consult Instructions for Use	-10°C	Temperature Limit (Storage)
	Reusable Do not throw away	0%	Humidity limitation (Operating and Storage
IPX2	Protected against vertically falling water drops when enclosure tilted up to 15°	SN	Serial number
EC REP	Authorized representative in the European Community	\sim	Date of manufacture
X	Waste electrical and electronic equipments must be disposed of in accordance with the local applicable regulations, not with domestic waste	C € 0123	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC. 0123 is Notified Body Number

Technical Specifications

Pulse Oximeter		
SpO2 Range	70% to 100%	
SpO2 Resolution	1%	
	90% to 100% range: ±2%;	
SpO2 Accuracy	70% to 89% range: ±3%	
	<70%: unspecified; complies with EN ISO80601-2-61	
Alarms	Low voltage indicator,power auto-off after be-be sound. Dual wavelength LED	
Method		
Pulae Rate Range	30 to 250 bpm	
Pulse Rate Resolution	1 bpm	
Pulse Accuracy	±3 bpm	
LED Wavelengths	Red: approximately 660nm; Infrared: approximately 905nm	
Optical output power	Less than 15mW	
Power Supply Requirements	Note: The Oximeter does not include batteries.	
Batteries	1.5V (AAA) alkaline batteryX2 (IEC Type LR03)	
Adaptable Voltage Range	2.6V~3.6V	
Pulse Oximeter		
Only SpO2 function works	±Less than 20mA	
Display Parameters	SpO2, Pulse Rate, Pulse Waveform Display, Pulae Amplitude Indicator and Low Voltage Indicator	

Environmental Operating Temperature 5 to 40°C (41 to 104°F) Limit Storage Tempeerature -10 to 40°C (14 to 104°F) Limit Atmospheric Pressure limitation (Storage, 86KPa~106KPa Transportation and Operating) Classification Medical device: Class II Equipment (II a EU Directive 93/42/EEC) Protection Against Liquids: IPX2 Dimension and Weight: 31.5g (Not including batteries) Weighting Size: 61*34*30.5mm Compliance Compliant with Item Safety Standards: EN 60601-1,EN60601-1-2 Equipment classification Type of protection nternally powered equipment (on battery power) Degree of protection Type BF Applied part Mode of operation Continuous Front panel and ISO15223-1 case labeling EN ISO80601-2-61 Pulse oximeter

Product parts and accessories

The Pulse Oximeter is consists display screen, power button, applied part, battery cover and sling. Detail of the instrument see Figure2.



• Do not have the oximeter immerged in liquid. When it needs cleaning, please wipe its surface with disinfect solution by soft material. Do not spray any liquid on the device directly

• When cleaning the device with water, the temperature should be less than 60 °C.

Unpacking and Inspection

Remove the instrument of PULSE OXIMETER from the shipping carton and examine for signs of shipping damage. Please check all materials against the packing list. Save the invoice, bill of lading and all packing materials. These may be required if it is necessary to process a claim with the carrier. If anything is missing or damaged, please contact the Technical Service Department. You can contact by:

- Phone: +86 755 61120085
- Fax: +86 755 61120055
- Email: user07@med-linket.com

Included in the package:

Description	Qty
PULSE OXIMETER (instrument)	1 Piece
User Manual	1 Piece
Sling	1 Piece

Section 2- Operation

Installation and Verification

 Battery installation Caution: The Pulse Oximeter does not operate with dead batteries and can not be powered by external power source does not input outer power. Install new batteries. Install new batteries.

- 1. Open and remove the battery cover.See Figure 3.
- 2. Insert two "AAA" size batteries, making sure the battery's positive and negative poles are correctly oriented in the battery compartment,
- 3. Install the battery cover. Battery cover



Figure 3: Installing Batteries

- WARNING: Explosion hazard. Do not use the PULSE OXIMETER in the ∕!∖ presence of flammable anesthetics mixed with air, with oxygen, or nitrous
 - oxide. WARNING: To ensure accurate performance and prevent device failure, do not expose the PULSE OXIMETER to extreme moisture such as rain.

6

Performance Verification 1.Performance Tests

- The power-up performance test verifies that the PULSE OXIMETER is ready for patient monitoring
- 2.Power-On Self-Test
- Before using the PULSE OXIMETER, you must verify that the PULSE OXIMETER is working properly and is safe to use. Proper working conditions are verified each time when the PULSE OXIMETER is turned on as described in the following procedure. The verification procedure (POST) takes 2 to 3 seconds to complete

Caution: If any indicator or display element does not light when the PULSE OXIMETER is turned on, do not use the PULSE OXIMETER. Instead, contact qualified service personnel, your local MED-LINKET representative, or MED-LINKET's Technical Services Department.

Note: Physiological conditions, medical procedures, or external agents that may interfere with the PULSE OXIMETER's ability to detect and show measurements, including dysfunctional hemoglobin, arterial dyes, low perfusion, dark pigment, and externally applied coloring agents such as nail polish, dye, or pigmented cream.

Note: The Pulse Oximeter automatically starts the Power-On Self-Test (POST) to ensure that its internal circuits are functioning properly.

Procedure

- Turn on the PULSE OXIMETER by pressing the Menu button.
- 2. After the device completes the Power-On Self-Test (POST), it will directly switch to measure interface.

General Operation

- Preparative for operating
 - 1) Open the battery cover carefully and then install two "AAA" Alkaline batteries according to the (+/-) electrodes.
- 2) Press the power button for 1 second to activate the device. SpO2 measure
 - 1) Open the clip of PULSE OXIMETER, See figure 4. ①.
 - 2) Place a finger (The preferred application is middle finger and index finger)





5.Get the information of SpO2 directly from screen display.



- The waste of PULSE OXIMETER must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the
- decommissioning of your equipment. Please refer to the correlative literature about the clinical restrictions and caution.
- This device is not intended for treatment.
- The LCD panel contains toxic chemicals. Do not ingest chemicals from a broken LCD panel.
- Do not modify this equipment without authorization of the manufacturer.

Latex Content Statement

About This Manual

understanding of operation.

Intended Use

measurement.

Section 1- Overview

About the Pulse Oximeter

Identification of Front Panel and Symbols

buttons, symbols, controls, displays and indicators.

98

connectors, switches.

does not exceed 35°C.

The PULSE OXIMETER and accessories are not made with natural rubber latex in any location that may result in patient contact.

1

The PULSE OXIMETER is to be operated by qualified personnel only. Before

The Pulse Oximeter is intended for continuous use or spot checking in measuring

in hospitals, physician's office, clinical settings and home care environment.

and displaying functional arterial oxygen saturation (SpO2)and pulse rate of patients

The device contains a dual light source (red LED and infrared red LED) and a photo

detector. Bone, tissue, pigmentation and venous vessels normally absorb a constant

amount of light over time. The arteriolar bed normally pulsates and absorbs variable

amounts of light during the pulsations. The ratio of light absorbed is translated in an

oxygen saturation measurement (SpO2). Because a measurement of SpO2 is

dependent on light from the device, excessive ambient light can interfere with this

Refer to the PULSE OXIMETER Operator's manual for a complete description of all

servicing this product, read the operator's manual carefully and a thorough

on the silicone (Ensure the finger position is correct that the LED $\,$ (irradiancy) window against finger prominence and the accepting window against finger lunula), see figure 4.2, and then clip the finger, see figure 4.33) Turn on the PULSE OXIMETER by pressing the Power button " (G) ". 4) Get the information of SpO2 directly from screen display.



Figure 4: measurement Note: When put finger into the silicone cushions of the clip, make sure nail is upturned.

7

Principle of Measurement

The measurement of PULSE OXIMETER is uses a multi-functional oxyhemoglobinometer to transmit some narrow spectrum light bands through blood samples, and to measure attenuation of spectrum with different wavelengths according to the characteristic that RHb, O2Hb, Met Hb and COHb absorb the light of different wavelength, thereby determining O2Hb saturation of different fractions. O2Hb saturation is called "fractional" O2Hb saturation

Oppositely, pulse oxygen oximeter measure functional O2Hb saturation:

Present Pulse oximeter transmits light of two wavelengths only, red light and infrared, to differentiate HbO2 from HbR. One side of the sensor contains two LEDs, and the other side contains a photoelectric detector. Pulse oximeter measures HbO2 saturation in the blood by the light plethysmograph when the pulse beats. The result is quite precise when HbO2 saturatiion is between 70% to 100%.

Clinical Restrictions

- 1) As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of the testee is required. For a testee with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO2 waveform will decrease. In this case, the measurement will be more sensitive to interference.
- 2) For those with a substantial amount of staining dilution drug such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (MetHb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO2 determination by this monitor may be inaccurate.
- 3) The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO2 measure.
- 4) The SpO2 value serves only as a reference value for judgement of anemic anoxia and toxic anoxia, some patients with serious anemia may also report good SpO2 measurement.

Attentions

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture
- If the oximeter gets wet, please stop using.
- When it is carried from cold environment to warm and humid environment, please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the oximeter is not permitted. Refer to User's Manual for instructions of cleaning and disinfection.

5

Safety 1) Safety

Instructions for safe operations

- Check the Pulse Oximeter periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected once a week at least. Please stop using the monitor when there is obvious damage.
- Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves
- At elevated ambient temperatures, patient skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41 °C on the skin if the initial skin temperature does not exceed 35°C
- Please remove the finger from the instrument to stop measure then the Pulse Oximeter will power off automatically within 8 seconds if the instrument must be closed for the urgent status.

2) Warnings <u>/</u>

- Explosive hazard—DO NOT use the Pulse oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- DO NOT use the oximeter while the testee is under measurement of MRI and CT.
- Be cautious of the sling. Please do not break the sling during usage to avoid device damage. Please don't use sling if allergic to it.
- Please don't use this product if you allergic to silicone cushion and ABS plastic.
- Please dispose the device, batteries and package materials (including plastic bag, foam and carton) according to local law.

3) The attention of Operation

- The equipment should be fully tested to see if it can be used normally before using.
- The finger should be placed properly (see figure 4 of this manual), otherwise it may cause inaccurate measurement.
- The SpO2 sensor and photoelectric receiving tube should be arranged in a way with the testee's arteriole in a position in between
- The SpO2 sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff of receiving intravenous injection
- Make sure the optical path is free from any optical obstacles like rubberized fabric; otherwise it may result in venous pulsation and inaccurate measure of SpO2.
- Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- Strenuous action of the testee or extreme electrosurgical interference may also affect the accuracy.
- Testee cannot use enamel or other makeup.
- Please clean and disinfect the device after operating according to the user manual.

Figure 1: PULSE OXIMETER Front Panel

78 0

1— Menu button/Power button	5— Waveform Display
2— %SpO2 Display	6— Bar graph (The Pulse Amplitude Indicator)
3— Low Battery indicator	7— Screen turn switch
4— Pulse Rate Display (bpm)	

8

Pulse Oximeter BZ1321B REV:A.1

Function Setting Introduction

Press the Pulse Oximeter Menu button to power on and access to the testing interface, or press the Menu button repeatedly during normal operation sequentially switch parameter-setting interfaces to set up the parameters and then return to the POST display. Settable parameters include high and low SpO2 limit, high and low bpm limits, high and pulse beep volume. The device will power off automatically within 8 seconds when there is no any signals input, andcan also use the menu button under parameter-setting interfaces to turn the PULSE OXIMETER off.

Menu Setting

Table 1: Instruction for Menu setting equipment and Figures Function Instruction for operation • Power on Turn on the PULSE OXIMETER by pressing the Power "on" and Menu/Power button ' G "off" • Power off setting Short press the button, move the cursor to © 98 78 0 select the item of "power off", and then long press the button to turn the power off. The device will power off automatically within 8 seconds when there is no any signal input. Setting enter Long press the button to enter the interface of settings. The setting interface of PULSE OXIMETER includes "Aim Setup 1", "Aim Setup 2" and "Sounds Setup". Setting enter and exit Alm Settings Alm Setupl on Beep off Power off Restore * OK Exit PULSE OXIMETER setting interface E xit • Short press the button, move the cursor to select the item of "Exit", long press the button return to the POST display. "Alm" on or off setting Short time presses the menu button to enter the interface of settings of "Alm Setup 1". Move the cursor select the item of "Alm", and then long press the button turn the functions on or off. " Alm " Alm Settings Alm Setupl Alm on Beever off Restore on or off setting Short press the button, move the cursor to select the item of "Exit", and then long press the button return to the POST display. E xit "Beep" on "Beep" on or off setting Alm Settings Alm Setupl or off Short press the button, move the cursor to on setting select the item of "Beep", and then long press $\begin{cases} Beep \\ OK \end{cases}$ $K = 0 \\ K \end{cases}$ the button to turn the functions on or off. Exit Default setting Short press the button, move the cursor to select the item of "Restore", then long press the button to returns the PULSE OXIMETER to factory default setting. Default setting After completing the setting, the interface will indicate "OK". Move the cursor to select the item of "Exit" by short press the button, and then long press the button to return to the POST display. 9 Alm Settings Alm Setupl Beep Power off Restore Settings Alm Setup1 Beep Power off Restore on on off * ^{0K} * E xit E xit SpO₂ High Limit setting Long press the button to enter the interface of settings of "Alm Setup 2". Short press the button, move the cursor to select the item of "SpO₂ Alm Hi", long press the button to adjust the parameter of SpO₂ in the scope of 50% to 100%. SpO₂ High Alm Settings Sp02 Alm Hi 100 Sp02 Alm Lo * 94 PR Alm Lo 50 PR Am Lo 50 Limit setting E xit \mbox{SpO}_2 Low Limit setting Short press the button in the interface of "Alm SpO₂ Low Alm Setting Sp02 Alm Hi Sp02 Alm Lo PR Alm Hi PR Alm Lo Limit * ¹⁰⁰ * ⁹⁴ ¹³⁰ ⁵⁰ Setup 2", move the cursor to select the item of "SpO₂ Alm Lo", long press the button to adjust the parameter of SpO₂ in the scope of 50% to setting E xit 100% Pulse Rate Short press the button in the interface of "Alm Alm Setup Sp02 Alm H Sp02 Alm L PR Alm Hi PR Alm Lo 100 94 * 130 * 50 (PR) High Setup 2", move the cursor to select the item of Limit "PR Alm Hi", long press the button to adjust the setting parameter of BPM in the scope of 5-250bpm. Pulse Rate Short press the button in the interface of "Alm 100 94 130 * Alm Setup SpO2 Alm H SpO2 Alm L PR Alm Hi PR Alm Lo Setup 2", move the cursor to select the item of Age (cm) (kg) Skin Tone (PR) Gender "PR Alm Lo", long press the button to adjust the Low Limit M 31 160 70 Light setting parameter of BPM in the scope of 5-250bpm. 165 55 2# M 24 3# F 22 160 45 After completed above setting, press the Return to 4# M 29 175 60 button switch to any interface of setting, move the POST © 👯 78 0 1 5# F 22 160 49 the cursor to select the item of "Exit" to return display 6# 19 160 45 to the POST display. F 21 162 54 8# M 34 192 102 Light (White) 9# F 27 178 58

A Warning

especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 4 hours.

Electromagnetic Environment

The PULSE OXIMETER is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the PULSE OXIMETER can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PULSE OXIMETER as recommended below, according to the maximum output of the communications equipment.

Warning:

- PULSE OXIMETER should not be used adjacent to or stacked with other
 - that if adjacent or stacked use is necessary, PULSE OXIMETER should be observed to verify normal operation in the configuration in which it will be

lse Ovimeter i	is intended for us	e in the electron	hagnetic environment specified below
e customer or vironment.	the user of the P	ulse Oximeter sho	buld assure thatitis used in such an
Emissions test	complia	nce Ele	ctromagnetic environment-guidance
emissions	Group 1	The P	ulse Oximeter uses RF energy for its
PR 11		emiss	ions are very low and are not likely to
		cause any interference in nearby electronic	
omissions	Class B	equip The P	ulse Oximeter is suitable for use in
PR11		dome	stic establishment and in
		voltag	e power supply network which
		suppli	es buildings used for domestic
	-		
Guidance	& Declaration -	Electromagnet	ic Immunity
T-1-1-4	0.11		
Table4-	-Guidance & De	claration — ele	ctromagnetic immunity
PULSE OXI low. The custo th an environr	METER is intended omer or the user of nent.	d for use in the ele f the PULSE Oxin	ectromagnetic environment specified neter should assure that it is used in
nunity test	IEC 60601	Compliance	Electromagnetic environment -
,	test ievel	ievel	guidance
rostatic	± 8 kV contact	± 8 kV contact	Floors should be wood, concrete or
arge (ESD)	$\pm 2 \text{ kV}, \pm 4 \text{ kV},$	$\pm 2 \text{ kV}, \pm 4 \text{ kV},$	ceramic tile. If floors are covered
1000-4-2	±8 KV,	± 8 KV, ±15KV ai	with synthetic material, the relative
	ттэ ку аlı.		numicity should be at least 30 %.
	101//m	10 \//m	
romagnotic	80MHz to 2 7CH-	10 1/11	
omagnetic	80% AM at 1kHz		
000-4-3			
			<u> </u>]
		12	
		12	
	20.4/m	12 20.4/m	Power frequency magnetic field-
er frequency	30 A/m 50 Hz or 60 Hz	12 30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields
r frequency tetic field 1000-4-8	30 A/m 50 Hz or 60 Hz	12 30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a twical location in a twical
r frequency etic field 1000-4-8	30 A/m 50 Hz or 60 Hz	12 30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital
r frequency etic field .000-4-8	30 A/m 50 Hz or 60 Hz	12 30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
r frequency etic field 1000-4-8	30 A/m 50 Hz or 60 Hz	12 30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
r frequency etic field 1000-4-8	30 A/m 50 Hz or 60 Hz	12 30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
r frequency letic field 1000-4-8 tion 5-	30 A/m 50 Hz or 60 Hz Measurem	12 30 A/m 50 Hz or 60 Hz ent Valida	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
r frequency letic field 1000-4-8 tion 5- 1	30 A/m 50 Hz or 60 Hz Measurem ter accuracy ha	12 ^{30 A/m} 50 Hz or 60 Hz ent Valida s been validate	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial
r frequency letic field 1000-4-8 tion 5- 1 ulse oxime sample ref	30 A/m 50 Hz or 60 Hz Measurem ter accuracy ha erence measure	12 30 A/m 50 Hz or 60 Hz ent Valida s been validate ed with a CO-O	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation
tion 5- i ulse oxime sample ref healthy ad	30 A/m 50 Hz or 60 Hz Measurem ter accuracy ha: erence measure ult volunteers w	12 ^{30 A/m} ^{50 Hz or 60 Hz} ent Valida s been validate ed with a CO-O ith saturation le	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation typels between 70% and 100% SaO2
er frequency hetic field 1000-4-8 tion 5- I ulse oxime sample ref healthy ad	30 A/m 50 Hz or 60 Hz Measurem ter accuracy has erence measure ult volunteers w	12 30 A/m 50 Hz or 60 Hz ent Validate ed with a CO-O ith saturation le	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2
tion 5- i ulse oxime sample ref healthy ad	30 A/m 50 Hz or 60 Hz Measurem ter accuracy hat erence measure ult volunteers w	12 30 A/m 50 Hz or 60 Hz ent Valida s been validate ed with a CO-O ith saturation le	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2
tion 5- i ulse oxime sample ref healthy ad studied.	30 A/m 50 Hz or 60 Hz Measurem ter accuracy has erence measure ult volunteers w	12 30 A/m 50 Hz or 60 Hz ent Valida s been validate ed with a CO-O ith saturation le	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2
frequency etic field 000-4-8 ion 5- ulse oxime ample ref healthy ad tudied. ct Demo	30 A/m 50 Hz or 60 Hz Measurem ter accuracy have erence measure ult volunteers w ographics haracteristics for	30 A/m 50 Hz or 60 Hz ent Valida s been validate ed with a CO-O ith saturation le	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2
r frequency etic field 1000-4-8 tion 5- 1 ulse oxime sample ref healthy ad tudied. ect Demo	30 A/m 50 Hz or 60 Hz Measurem ter accuracy have erence measure ult volunteers w ographics haracteristics for	12 30 A/m 50 Hz or 60 Hz ent Valida s been validate ed with a CO-O ith saturation le	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2 as follow table 5
r frequency etic field .000-4-8 tion 5- 1 ulse oxime sample ref healthy ad tudied. ct Demo pulation c pulation c	30 A/m 50 Hz or 60 Hz Measurem ter accuracy hater erence measure ult volunteers w ographics haracteristics fo LSE OXIMETEF	12 30 A/m 50 Hz or 60 Hz ent Valida ent Validate ed with a CO-O ith saturation le r those studies & Clinical study	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2 as follow table 5 Subject Demographics Record.
ion 5- A lise oxime sample ref healthy ad tudied. ct Democ pulation c ble 5-PUl	30 A/m 50 Hz or 60 Hz Measurem ter accuracy hater erence measure ult volunteers w ographics haracteristics for LSE OXIMETER	12 30 A/m 50 Hz or 60 Hz ent Validate ed with a CO-O ith saturation le r those studies R Clinical study abt Weight	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. tion d in human studies against arterial ximeter. In a controlled desaturation evels between 70% and 100% SaO2 as follow table 5 Subject Demographics Record.

Light

Medium Dark

Light (White)

Light (White)

Dark dark

Dark dark

Light

Light

Light

78

The final analysis was performed on 249 data points collected across 11 subjects.

24 174 80 Dark dark

23 178

26 169 65

10#

11#

12#

M

F

M

ARMS Results:

Bland-Altman Plot Comparing the SpO2 Difference between the PULSE b) OXIMETER and the Reference CO-Oximeter During Non-Motion Conditions



ection 6- Service and Maintenance

leaning and Disinfecting

- Clean the surface of the oxmeter by using a soft cloth dampened with either a commercial, non-abrasive cleaner or a solution of 70% isopropyl alcohol in water, and wiping it lightly the surfaces of the oximeter.
- Please switch off pulse oximeter before cleaning. Clean the LED and Please switch oil puise oxinities of sole of the photo-sensor with moist cloth or cotton ball and alcohol gently.
- The aforementioned general cleaning process is not for infection prevention. Please contact the specialist for process of contagious infection.

alibrating

Please use the SpO2 simulator of Fluke Biomedical index 2 to calibrate PULSE OXIMETER for the function of SpO2 measure. The calibration must be operated to by qualified personnel only. The SpO2 accuracy can be validated in human studies against arterial blood sample reference measured with a CO-oximeter. All of the process of the clinical study must becomplied with standard of EN ISO80601-2-61.

epairing and Maintenance

1.	Please change the batteries when the low-voltage indicator lightens.
2.	Please clean the surface of the device before using. Wipe the device with
	alcohol first, and then let it dry in air or clean it by dry clean fabric.

Please take out the batteries if the oximeter is not in use for a long time.

15

• Uncomfortable or painful feeling may appear if use the device ceaselessly,

• For the special patients, there should be a more prudent inspecting in the placing process. The device cannot be clipped on the edema and ter

The SpO2 accuracy performance of each pulse oximeter and sensor comb

Asian (Chinese)

Asian (Chinese)

Asian (Chinese)

Asian (Chinese)

Asian (Chinese)

Asian (Chinese)

Caucasian

Caucasian

Caucasian

African

African

African

- The best storage environment of the device is -10°C to 40°C ambient temperature and not higher than 80% relative humidity.
- Please maintain properly for ensuring the device can be used normally. The device needs to be calibrated once a year (or according to the
- calibrating program of hospital). It can also be performed at state-appointed agent or just contact us for calibration.

/ Warnings

- High-pressure sterilization cannot be used on the device.
- Do not immerse the device in liquid.
- It is recommended that the device should be kept in a dry environment. Humidity may reduce the using life , or even damage the device.

isposal

- Used batteries should not be disposed of in the household rubbish. Used Batteries should be deposited at a collection point.
- At the end of its life, the appliance should not be disposed of in household rubbish. Enquire about the options for environment-friendly and appropriate disposal. Take local regulations into account.

Warranty

Our company warrants pulse oximeter at the time of its original purchase and for the subsequence time period of one year. The warranty does not cover the followings:

17

- The device series number label is torn off or cannot be recognized.
- Damage to the device resulting from misconnection with other devices. ۲
- Damage to the device resulting from accidents. Changes performed by users without the prior written authorization of the ۲
- company.

Qualified certificate

 The light (the infrared is invisible) emitted from the device is narmful to the eyes, so the user and the maintenance man should not stare at the light. Testee cannot use enamel or other makeup. Testee's fingernail cannot be too long. Please refer to the correlative literature about the clinical restrictions and caution. This device is not intended for treatment. The user is not allowed to repair the equipment. Changes or modification not expressly approved by Shenzhen Med-link may void the warranty. Removing the batteries to avoid battery leakage and device damage if long time no use. 	Arms = $\sqrt{\sum_{i=1}^{n} (SpO_{2i} - RefSaO_{2i})^2}$ Where: ARMS is the accuracy root mean square. SpO2 is the test pulse oximeter readings during sample i. RefSaO2 is the Average Reference CO-Oximeter functional oxygen saturation reading during sample i. n is the number of points. The detail of the ARMS Results is below table 6 and table 7.	PRODUCT NAME See product labels PRODUCT MODEL See product labels DATE INSPECTOR	
10	13	16	
<section-header><section-header><section-header><section-header><section-header><section-header><section-header><section-header><section-header><section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header></section-header>	Table 6—overall Average Root Mean Square (ARMS) for PULSE OXIMETER in the SpO2 range of 70%-100%. Image: Compared to Avg. Reference Origination of the compared to Avg. Reference PULSE OXIMETER in a clinical study. Specification 70-100% Areas Pulse Oximeter Image: Compared to Avg. Reference Co-Oximeter, Functional SaO2 Apr 6-8, 2012 SaO2 ranges of 80-90% Areas Pulse Oximeter SaO2 ranges of 80-90% Areas Pulse Oximeter SaO2 ranges of 90-100% Areas Pulse Oxim	Shenzhen Med - link Electronics Tech Co., Ltd. 4th and 5th Floor, Building Tow,Hualian Industrial Zone, Xinshi Community,Dalang Stree,longhua 518109 Shenzhen,PEOPLE'S REPUBLIC OF CHINA Email: sales@med - linket.com Shanghai International Holding Corp. GmbH (Europe) . Ad: Eiffestrasse 80, 20537 Hamburg, Germany Tel: +49-40-2513175 Fax: +49-40-255726	
This device has been tested and found to comply with the limits for medical devices to the IEC/CN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare environments (for example, electrosurgical units, cellular phones, mobile two-way radios, electrical appliances, and high-definition television), it is possible that high levels of such interference due to close proximity or strength of a source may result in disruption of performance of this device. This Fingertip pulse oximeter is not designed for use in environments in which the pulse can be obscured by electromagnetic interference. During such interference, measurements may seem inappropriate or the monitor may not seem to operate correctly.	65 70 75 80 85 90 95 100 Item 70100 90100 80<90		