


New Product



Temperature Probe

Product Advantage

- ★ Made with high-quality temperature sensor, accurate measurement (25~45°C, accuracy $\pm 0.1^{\circ}\text{C}$), quick response;
- ★ It has a wide range of applications and can be used for adult, pediatric and neonatal patients;
- ★ Latex free, good biocompatibility, no allergy to human body.

Scope of Application

It can be used with a multi-parameter monitor to transmit the patient's temperature signal.

Product Parameter

Compatible Brand	Biolight A series (A2/A3/A5/A6/A8), Q series (Q3/Q4/Q5/Q6/Q7)		
photo	OEM#	Order Code	Specification
A	15-031-0005	W0113A	2pin, Adult Skin-surface, 10ft
B	15-031-0012	W0113B	2pin, Adult Rectal/Esophageal, 10ft
C	/	W0113C	2pin, Pediatric Skin-surface, 10ft
D	/	W0113D	2pin, Pediatric Rectal/Esophageal, 10ft
E	/	W0113G	Temperature Adapter Cable, 8.2ft (2.5m)
F	8001642	W0001E	Disposable Temperature Probe Skin-surface, 32in, 24pcs/box

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address: 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com