


Medlinket Products



Direct-Connect ECG Cables

Product Name: Direct-Connect ECG Cables

Brand: Medlinket

Patient end terminal: Snap, Clip, Grabber

Applicable Person: Adult, Pediatric

OEM/ODM: Yes

Quality Certification: CE, FDA, NMPA

*For more product details, check out the information below or contact us directly, marketing@med-linket.com

Compatibility Information:

Compatible Brand	Original Model
Biolight	15-031-0014, 15-031-0004
Comen	040-000320-00(AHA), 040-000368-00(AHA), 040-000321-00(AHA), 040-000369-00(AHA)
Drager/Siemens	/
Datex	/
GE-Marquette>Hellige	/
Mindray	/
Nihon Kohden	/
Philips	/
Schiller	/
Spacelabs	/

Specification:

- 1) Leads: 3LD, 4LD, 5LD, 6LD
- 2) Standard: AHA, IEC



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

- 3) Patient end terminal: Snap, Clip, Grabber
- 4) Material: TPU, PVC
- 5) Package: 1pcs/bag
- 6) Defibrillation with stand: 1K Ω

Main Instrument Connector:



Patient End Terminal:



Product Advantages:

- 1. For adult and pediatric patients;
- 2. Meeting the requirements of EC53;
- 3. Integrated design, convenient for use and maintenance;
- 4. Flexible and durable cables, enduring repeated cleaning and disinfection;
- 5. Outstanding shielding property and anti-interference performance, protecting ECG signal from being interfered;
- 6. The electrode connector has a clear color mark, which can quickly and easily connect the electrode position;
- 7. Latex free.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of SpO₂, temperature, EEG, ECG, blood pressure, EtCO₂, high-frequency electrosurgical products, etc. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Direct-Connect ECG Cables, ECG Cables, manufacturer, factory

Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com