




Disposable Defibrillation Pads

Product Advantages

- 1.Integrates the functions of defibrillation, pacing, monitoring, and cardioversion.
- 2.The pre-connected cable saves time and helps quickly start the right therapy , such as CPR or defibrillation.
- 3. The pads with color - coded operation indicators help to rapidly position and simplify the operation.
- 4. Low contact impedance to reduce skin burns.

Suggested Applications

Defibrillation, pacing, monitoring, and cardioversion

Declaration: Shenzhen Med-Link Electronics Tech Co., Ltd. ("MedLinket") hereby declares that all registered trademarks, product names, model numbers, etc., mentioned in this document belong to respective original owners or manufacturers.The information referenced herein is intended to clearly describe the compatibility features of MedLinket products for reference only and should not be construed as an operational guideline for medical institutions or related entities. If any adverse consequences cause from misunderstandings or incorrect use of the information, MedLinket shall not be held liable.

Ordering Information

Compatible Models	Connector Pictures	OEM #	Order Code	Description
CU Medical System Paramedic CU-ER1 Paramedic CU-ER2 Paramedic CU-ER3 Paramedic CU-ER5 con cavo adattatore ER-5 OA03 Paramedic CU-ERT			DE0001A	Adult/Pediatric≥25 kg, 1.2m
			DE0001P	Pediatric or Infant <25kg, 1.2m
Medtronic Physio-Control LIFEPAK 10/12/20/500/1000 Osatu Bexen REANIBEX-200/300/500/700/800 Cardioline ELIFE 700 Mindray BeneHeart D1/D2/D3/D5/D6 Comen F1/F1A/F2/F2A		0651-30-77007 (Mindray)	DE0002A	Adult/Pediatric≥25 kg, 1.2m
			DE0002P	Pediatric or Infant <25kg, 1.2m
CMICS AED-201 RadianQbio HR-501/503/501T/ 503T/701T/ 701T Plus			DE0003A	Adult/Pediatric≥25 kg, 1.2m
			DE0003P	Pediatric or Infant <25kg, 1.2m
CMOS DRAKE Dea Life 400 Futura(after 2018) Cardioversor Vivo Gold (after 2024)		F7988/CM (CMOS DRAKE)	DE0004A	Adult/Pediatric≥25 kg, 1.2m
			DE0004P	Pediatric or Infant <25kg, 1.2m
CU Medical System IPAD NF1200/1201 Ambulanc Amoul i3 Amoul i5 CMOS DRAKE Dea Life 400 Futura (before 2018) Cardioversor Vivo Gold (before 2024)		F7959/W (CMOS DRAKE)	DE0005A	Adult/Pediatric≥25 kg, 1.2m
			DE0005P	Pediatric or Infant <25kg, 1.2m
Zoll PD 1200/1400/1600/1700/2000 E Series M Series R Series M&B AED7000			DE0006A	Adult/Pediatric≥25 kg, 1.2m
			DE0006P	Pediatric or Infant <25kg, 1.2m
Philips 940010XX9400940020XX & 9420XX &9400 M1722A/B M1723A/B M1724A M2475B/E/S/EM Series			DE0007A	Adult/Pediatric≥25 kg, 1.2m
			DE0007P	Pediatric or Infant <25kg, 1.2m
Philips HeartStart FR2/FR2+ HeartStart FR3 HeartStart FRx HeartStart XL+; M3508A plug connector cable HeartStart MRx; M3508A plug connector cable HeartStart Intrepid; M3508A plug connector cable Efficia DFM100; M3508A plug connector cable		M3713A	DE0008A	Adult/Pediatric≥25 kg, 1.2m
		M3717A, M3870A	DE0008P	Pediatric or Infant <25kg, 1.2m
Nihon Kohden AED-21XX series		H324B	DE0009A	Adult/Pediatric≥25 kg, 1.2m
			DE0009P	Pediatric or Infant <25kg, 1.2m
Primedic HeartSave PAD HeartSave AED HeartSave AED-M HeartSave 6 HeartSave 6S		96389	DE0010A	Adult/Pediatric≥25 kg, 1.2m
			DE0010P	Pediatric or Infant <25kg, 1.2m



Shenzhen Med-link Electronics Tech Co.,Ltd

Address:Zone A of 1st and 2nd Floor, and 3rd Floor, Building A, No. 7, Tongsheng Industrial Park Road, Shanghengliang Community,

Dalang Street, Longhua District, 518109 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Tel: 86-400-058-0755 / 86-0755-61120299

Email: marketing@med-linket.com

Website: www.med-linket.com

