



Ordering Information

Product Images	Order Code	Dimensions	Description	Sterilization	Hypoallergenic	Radiolucent	Shelf Life	Population	Clinical Scenarios	Packaging
	V0015A	φ50	4mm Carbon Snap , Cloth	×	×	✓	2 Years	Adult	CT(X-ray) DR(X-ray) DSA(X-ray) MRI(Top Recommended)	25pcs/bag, 250pcs/box
	V0014A	φ50	4mm Metal Snap , Cloth	×	×	×	2 Years	Adult	Suitable for hospital - wide use, especially in the ER, OR, Catheter Lab, ICU and for use by Paramedics.	25pcs/bag, 250pcs/box
	V0014NP	Φ40		×	×	×	2 Years	Pediatric		25pcs/bag, 250pcs/box
	V0014C	Φ30		×	×	×	2 Years	Pediatric Infant Neonate		25pcs/bag, 250pcs/box
	V0014AP	φ50	4mm Metal Snap , Foam	×	×	×	2 Years	Adult		25pcs/bag, 250pcs/box
	V0014CP (V0014C-I)	φ25		×	×	×	2 Years	Pediatric Infant Neonate		25pcs/bag, 250pcs/box
	V0014DP (V0014D)	10*25		×	×	×	2 Years	Infant Neonate		10pcs/bag, 500pcs/box

Declaration: Shenzhen Med-Link Electronics Tech Co., Ltd. ("MedLinket") hereby declares that all registered trademarks, product names, model numbers, etc., mentioned in this document belong to respective original owners or manufacturers.The information referenced herein is intended to clearly describe the compatibility features of MedLinket products for reference only and should not be construed as an operational guideline for medical institutions or related entities.If any adverse consequences cause from misunderstandings or incorrect use of the information, MedLinket shall not be held liable.

Disposable ECG Electrodes (Standard Type)

Product Features

- 1.Recommended for patients of all ages, including neonates, pediatrics, and adults.
- 2.Applicable in diverse medical scenarios, including diagnosis, monitoring, CT, DR, DSA, and MRI.
- 3.High-quality medical pressure-sensitive adhesives are difficult to peel off even if water or medical solutions contact the electrodes.
- 4.Adopts unique polymerization technology to reduce skin irritation for patients with sensitive skin.
- 5.Not manufactured with natural rubber latex, plasticizers, or mercury.

Offers OEM services to meet diverse client requirements, including electrodes of different materials, shapes, and patterns.



Shenzhen Med-link Electronics Tech Co.,Ltd
Address:Zone A of 1st and 2nd Floor, and 3rd Floor, Building A, No. 7, Tongsheng Industrial Park Road, Shanghenglang Community, Dalang Street, Longhua District, 518109 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
Tel: 86-400-058-0755 / 86-0755-61120299
Website: www.med-linket.com

