


Medlinket Products



Disposable ECG Leadwires

Product Name: Disposable ECG Leadwires

Brand: Medlinket

Applicable Person: Adult, Pediatric and neonate patients

OEM/ODM: Yes

Quality Certification: CE, FDA, NMPA

*For more product details, check out the information below or contact us directly, marketing@med-linket.com

Compatibility Information:

Compatible Brand	Original Model
Covidien	33103, 33105, 33105E, 33111, 33136R36
DIN Type Connector	M3915A (PHILIPS), 900716-001 (GE)
Drager	MS14556, MS14555, MP00877, MP00875, MS14560, MS14559, MP00881, MP00879, MS14682, MS14683, MP03123, MP03122
Datex	/
GE	E9008LF, E9008LH, E9003CL, E9003CN, E9008KB, E9008KD, E9002ZW, E9002ZZ
Mindray	0010-30-42734, 0010-30-42733, 0010-30-42731, 0010-30-42732, 0010-30-42735, 0010-30-42736, 0010-30-42727, 0010-30-42730
Philips	M1673A, M1674A, 989803173121, 989803174201, M1644A, M1645A, 989803173131, 989803174211, M1604A, M1602A, M1978A, M1976A
Siemens	/
Spacelabs	/

Specification:


- 1) Leads: 3LD, 5LD, 6LD
- 2) Standard: AHA, IEC
- 3) Patient end terminal: Snap, Grabber



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
 Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
 Email: sales@med-linket.com | Website: www.med-linket.com

Main Instrument Connector:



Patient End Terminal:



Product Advantages:

1. Single patient use: decreases the risk of cross contamination and improves patient care;
2. Fully shielded seadwires design: reduces the risk of Electromagnetic Interference (EMI);
3. Peelable ribbon cable design: prevents lead wire entanglement and can be customized to fit any patient's body size;
4. Side button and visual connection design: (1) Provide clinicians with a locking and visualization mechanism to achieve faster, efficient and firm connections; (2) Clinically proven to reduce the risk of "leads off" false alarms;
5. Easy to used electrode colors Lightweight smooth design: (1) Easy and quick lead placement; (2) Increase patient comfort.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of SpO₂, temperature, EEG, ECG, blood pressure, EtCO₂, high-frequency electrosurgical products, etc. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Disposable ECG Leadwires, ECG Cables, manufacturer, factory

Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com