


# New Product



## Disposable Suction Tube Accessories

### Product Advantage

- ★ Polymer material, anti-negative pressure, good flexibility;
- ★ Standard connector range, can fit a variety of models to facilitate the clinical use;
- ★ Non-sterile provision, disposable use, avoiding cross-infection;
- ★ Latex-free, cost-effective .

### Scope of Application

After matching with suitable equipment, it is used for drainage and suction of blood and waste fluid during and after the operation.

### Product Parameter

Order Code	Photo	ID Tube (mm)	Length(m)	Connector type	Package	Price Code
G5018A		5.0	1.8	Φ 7.7mm horn type	50pcs /box	H9/box
G5030A			3.0			J5/box
G5040A			4.0			K0/box
G5060A			6.0			M6/box
G5018B		5.0	1.8	Φ 7.7mm straight type	50pcs /box	H9/box
G5030B			3.0			J5/box
G5040B			4.0			K0/box
G5060B			6.0			M6/box



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
 Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
 Email: sales@med-linket.com | Website: www.med-linket.com

# New Product

Order Code	Photo	ID Tube (mm)	Length(m)	Connector Type	Package	Price Code
G5018AC		5.0	1.8	Φ7.7mm horn & elbow Type	50pcs /box	H9/box
G5030AC			3.0			J5/box
G5040AC			4.0			K0/box
G5060AC			6.0			M6/box
G5018BC		5.0	1.8	Φ7.7mm elbow & straight type	50pcs /box	H9/box
G5030BC			3.0			J5/box
G5040BC			4.0			K0/box
G5060BC			6.0			M6/box
G6518D		6.5	1.8	Φ9.9mm horn type	50pcs /box	J5/box
G6530D			3.0			L8/box
G6540D			4.0			O2/box
G6560D			6.0			S1/box
G6518E		6.5	1.8	Φ9.9mm straight type	50pcs /box	J5/box
G6530E			3.0			L8/box
G6540E			4.0			O2/box
G6560E			6.0			S1/box
G6518DF		6.5	1.8	Φ9.9mm horn & elbow Type	50pcs /box	J5/box
G6530DF			3.0			L8/box
G6540DF			4.0			O2/box
G6560DF			6.0			S1/box
G6518EF		6.5	1.8	Φ9.9mm elbow & straight type	50pcs /box	J5/box
G6530EF			3.0			L8/box
G6540EF			4.0			O2/box
G6560EF			6.0			S1/box



## Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
 Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community,  
 Dalang Street, Longhua District, 518109 Shenzhen, China  
 Email: sales@med-linket.com | Website: www.med-linket.com