







## 1 Disposable Temperature Probes

| Code Rule Description                                                                                                                                                                                                                                                                                        | Probe code                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Product Features                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>W <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/></p> <p>Probe type:<br/>E: Disposable Skin-Surface Probes<br/>F: Disposable Esophageal/Rectal Probes</p> <p>Compatible Brand code a 0001 indicated YSI 400 series</p> <p>Temperature Probes</p> | <p>E </p> <p>P </p> <p>F </p> <p>L </p> <p></p> <p></p> | <p>1. Single patient use, prevent cross infection.</p> <p>2. The accuracy of the thermistor is <math>\pm 0.1^{\circ}\text{C}</math> from <math>25^{\circ}\text{C}</math> to <math>45^{\circ}\text{C}</math>.</p> <p>3. Electromagnetic compatibility.</p> <p>4. Biocompatibility passed.</p> <p>5. The reflector can effectively insulate ambient temperature and radiant light (Skin-surface).</p> <p>6. The medical blue PVC tube provides atraumatic insertion. Leakage current of the temperature probe (thermistor and PVC tube) is well within clinically accepted non-hazardous levels (Esophageal/Rectal one).</p> <p>Disposable Skin-surface Probes</p> <p>Disposable Esophageal/Rectal Probes</p> <p>Disposable Esophageal/Rectal Probes</p> <p>Disposable Esophageal/Rectal Probes</p> <p>Disposable Ear Canal Probes canal probes</p> <p>Foley Catheter Temperature Probes</p> |

| OEM Manufacturer & Models                                                                                                                                              | Instrument Connector                                                                | Probe Type                                                                          | OEM #   | Order Code                                                         | Dimension | Descriptions                                                                                                  | Price |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------|--------------------------------------------------------------------|-----------|---------------------------------------------------------------------------------------------------------------|-------|
| Compatible YSI 400                                                                                                                                                     |                                                                                     |                                                                                     |         |                                                                    |           |                                                                                                               |       |
| Datex Ohmeda<br>Philips<br>Draeger<br>Edan<br>GE Marquette<br>MEK<br>Medtronic<br>Mindray ><br>Datascope<br>Nihon Kohden<br>Spacelabs<br>Welch Allyn<br>YSI 400 Series |    |    | 8001642 | W0001E                                                             | 28.8*30mm | 32in.(0.8m),Adult/Pediatric/Infant,Skin-surface Probes<br>24pcs/pk,25°C,<br>2.252KΩ                           |       |
|                                                                                                                                                                        |                                                                                     |    | 4499    | W0028E                                                             | φ30mm     | 32in.(0.8m),Adult/Pediatric/Infant,Skin-surface Probes<br>24pcs/pk,25°C,<br>2.252KΩ                           |       |
|                                                                                                                                                                        |                                                                                     |    | 4499    | W0001Q                                                             | 10*10mm   | 32in.(0.8m),Infant,Skin-surface,24pcs/pk,25°C,<br>2.252KΩ                                                     |       |
|                                                                                                                                                                        |                                                                                     |    | 8001644 | W0001P                                                             | 7FR       | 19in.(0.48m),Rectal/Esophageal 24pcs/pk,25°C,<br>2.252KΩ                                                      |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0001F                                                             | 9FR       |                                                                                                               |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0001L                                                             | 12FR      |                                                                                                               |       |
|                                                                                                                                                                        |                                                                                     |    | 8001644 | W0099P                                                             | 7FR       | 24in.(0.6m),Rectal/Esophageal 24pcs/pk,<br>25°C,2.252KΩ                                                       |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0099F                                                             | 9FR       |                                                                                                               |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0099L                                                             | 12FR      |                                                                                                               |       |
|                                                                                                                                                                        |                                                                                     |    | 8001644 | W0104P                                                             | 7FR       | 24in.(0.6m),Rectal/Esophageal 24pcs/pk,25°C,<br>2.252KΩ                                                       |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0104F                                                             | 9FR       |                                                                                                               |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0104L                                                             | 12FR      |                                                                                                               |       |
|                                                                                      |                                                                                     | W0101F-A                                                                            |         | 32in.(0.8m),Ear Canal24pcs/pk,<br>25°C,2.252KΩ                     |           |                                                                                                               |       |
|                                                                                     |                                                                                     | W0101F-P                                                                            |         | 32in. (0.8m), Ear Canal Probe 24pcs/pk, Pediatric,<br>25°C,2.252KΩ |           |                                                                                                               |       |
|                                                                                     |                                                                                     | W0109F                                                                              | 9FR     | 32in. (0.8m), Esophageal/Rectal 24pcs/pk, Adult 25°C,2.252KΩ       |           |                                                                                                               |       |
| Compatible YSI Series                                                                                                                                                  |                                                                                     |                                                                                     |         |                                                                    |           |                                                                                                               |       |
| YSI 400                                                                                                                                                                |  |  |         | W0103E                                                             | 28.8*30mm | 32in. (0.8m), 6.3mm mono plug<br>Skin-surface Probes<br>24pcs/pk,Adult/<br>Pediatric/Infant,<br>25°C, 2.252KΩ |       |
| Compatible YSI 10K Series                                                                                                                                              |                                                                                     |                                                                                     |         |                                                                    |           |                                                                                                               |       |
| Suitable for China<br>OEM patient monitor:<br>Tianrong:TR900D/E;<br>Anke: ASC553A3,<br>ASC553;Comen:<br>star 8000A/B/C, star<br>5000, 5000B/C                          |  |  |         | W0004E                                                             | 28.8*30mm | 32in.(0.8m),Skin-surface Probes 24pcs/pk,25°C,10KΩ<br>Adult/Pediatric/Infant                                  |       |
|                                                                                                                                                                        |                                                                                     |  |         | W0004P                                                             | 7FR       | 19in.(0.6m),<br>Rectal/Esophageal<br>24pcs/pk,25°C,10KΩ                                                       |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0004F                                                             | 9FR       |                                                                                                               |       |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0004L                                                             | 12FR      |                                                                                                               |       |

| OEM Manufacturer & Models                                                                                                                                       | Instrument Connector                                                                | Probe Type                                                                          | OEM #         | Order Code | Dimension | Descriptions                                                                                       | Price |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------|------------|-----------|----------------------------------------------------------------------------------------------------|-------|
| Compatible Mindray (China)                                                                                                                                      |                                                                                     |                                                                                     |               |            |           |                                                                                                    |       |
| PM6800 Series                                                                                                                                                   |    |    | 0011-30-37391 | W0013E     | 28.8*30mm | Skin-surface Probes<br>24pcs/pk, 25°C, 2.252KΩ                                                     |       |
| BENEVISION/T1/T5/T8<br>IMEC series<br>IPM Series/9800                                                                                                           |                                                                                     |    |               | W0013F     | 9FR       | 19in. (0.48m),<br>Esophageal/Rectal Probes<br>24pcs/pk, Adult/<br>Pediatric,<br>25°C, 2.252KΩ      |       |
| Compatible Philips                                                                                                                                              |                                                                                     |                                                                                     |               |            |           |                                                                                                    |       |
| All Philips temperature monitors that accept two-prong plugs                                                                                                    |    |    | 21091A        | W0003E     | 28.8*30mm | 32in.(0.8m),Skin-surface Probes 24pcs/pk,25°C, 2.252KΩ                                             |       |
|                                                                                                                                                                 |                                                                                     |    | 21090A        | W0003P     | 7FR       | 19in.(0.48m),Rectal/Esophageal 24pcs/pk,25°C, 2.252KΩ                                              |       |
|                                                                                                                                                                 |                                                                                     |                                                                                     |               | W0003F     | 9FR       |                                                                                                    |       |
|                                                                                                                                                                 |                                                                                     |                                                                                     |               | W0003L     | 12FR      |                                                                                                    |       |
| Compatible COMEN                                                                                                                                                |                                                                                     |                                                                                     |               |            |           |                                                                                                    |       |
| NC series<br>(NC3/5/8/10/12/19)<br>STARseries<br>(STAR5000A/B/C/D/E/F/H )<br>STAR8000<br>(STAR8000A/B/C/D/E/F/H)<br>C series<br>(C20/30/50/60/70/80/86/90 /100) |   |   |               | W0108F     | 9FR       | 24in. (0.6m),<br>Esophageal/RectalProbes<br>Adult/Pediatric,<br>24pcs/pk,<br>25°C, 10KΩ            |       |
| Compatible Med-linklet                                                                                                                                          |                                                                                     |                                                                                     |               |            |           |                                                                                                    |       |
| AM-807-C<br>AM-807B<br>AM-807B-C                                                                                                                                |  |  |               | W0024E     | 28.8*30mm | 32in. (0.8m),<br>Skin-surface Probes<br>24pcs/pk,Adult/Pediatric/Infant /Neonatal<br>25°C, 2.252KΩ |       |
|                                                                                                                                                                 |  |  |               | W0114E     | 28.8*30mm | 32in. (0.8m),<br>Skin-surface Probes<br>24pcs/pk,Adult/Pediatric/Infant /Neonatal<br>25°C, 2.252KΩ |       |
| AM-807<br>AM-808                                                                                                                                                |  |  |               | W0024F     | 9FR       | 19in. (0.48m),<br>Esophageal/Rectal Probes<br>24pcs/pk, Adult/<br>Pediatric,<br>25°C, 2.252KΩ      |       |

| OEM Manufacturer & Models                                                                                                                                              | Instrument Connector                                                                | Probe Type                                                                          | OEM #   | Order Code | Dimension | Descriptions                                                                                      | Pice |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------|------------|-----------|---------------------------------------------------------------------------------------------------|------|
| <b>Compatible YSI 400</b>                                                                                                                                              |                                                                                     |                                                                                     |         |            |           |                                                                                                   |      |
| Datex Ohmeda<br>Philips<br>Draeger<br>Edan<br>GE Marquette<br>MEK<br>Medtronic<br>Mindray ><br>Datascope<br>Nihon Kohden<br>Spacelabs<br>Welch Allyn<br>YSI 400 Series |    |    | 8001642 | W0001ES    | 28.8*30mm | Sterilized, 32in.(0.8m),Adult/Pediatric/Infant,Skin-surface Probes 24pcs/pk,25°C, 2.252KΩ         |      |
|                                                                                                                                                                        |                                                                                     |    | 4499    | W0028ES    | φ30mm     | Sterilized, 32in.(0.8m),Adult/Pediatric/Infant,Skin-surface Probes 24pcs/pk,25°C, 2.252KΩ         |      |
|                                                                                                                                                                        |                                                                                     |    | 4499    | W0001QS    |           | Sterilized, 32in.(0.8m),Infant,Skin-surface,24pcs/pk,25°C,2.252 KΩ                                |      |
|                                                                                                                                                                        |                                                                                     |    | 8001644 | W0001PS    | 7FR       | Sterilized, 19in.(0.48m),Rectal/Esophageal 24pcs/pk,25°C, 2.252 KΩ                                |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0001FS    | 9FR       |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0001LS    | 12FR      |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |    | 8001644 | W0099PS    | 7FR       | Sterilized, 24in.(0.6m),Rectal/Esophageal 24pcs/pk, 25°C,2.252KΩ                                  |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0099FS    | 9FR       |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0099LS    | 12FR      |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |    | 8001644 | W0104PS    | 7FR       | Sterilized, 24in.(0.6m),Rectal/Esophageal 24pcs/pk,25°C, 2.252KΩ                                  |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0104FS    | 9FR       |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0104LS    | 12FR      |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |   |         | W0101FS-A  |           | Sterilized, 32in.(0.8m),Ear Canal, 24pcs/pk, Adult, 25°C,2.252KΩ                                  |      |
|                                                                                                                                                                        |                                                                                     |  |         | W0101FS-P  |           | Sterilized, 32in. (0.8m), Ear Canal Probe 24pcs/pk, Pediatric, 25°C,2.252KΩ                       |      |
|                                                                                                                                                                        |                                                                                     |  |         | W0109FS    | 9FR       | Sterilized, 32in. (0.8m), Esophageal/Rectal 24pcs/pk, Adult 25°C,2.252KΩ                          |      |
| <b>Compatible YSI Series</b>                                                                                                                                           |                                                                                     |                                                                                     |         |            |           |                                                                                                   |      |
| YSI 400                                                                                                                                                                |  |  |         | W0103ES    | 28.8*30mm | 32in. (0.8m), 6.3mm mono plug Skin-surface Probes 24pcs/pk,Adult/ Pediatric/Infant, 25°C, 2.252KΩ |      |
| <b>Compatible YSI 10K Series</b>                                                                                                                                       |                                                                                     |                                                                                     |         |            |           |                                                                                                   |      |
| Suitable for China<br>OEM patient monitor:<br>Tianrong:TR900D/E;<br>Anke: ASC553A3,<br>ASC553;Comen:<br>star 8000A/B/C, star<br>5000, 5000B/C                          |  |  |         | W0004ES    | 28.8*30mm | 32in.(0.8m),Skin-surface Probes 24pcs/pk,25°C,10KΩ Adult/Pediatric/Infant                         |      |
|                                                                                                                                                                        |                                                                                     |  |         | W0004PS    | 7FR       | 19in.(0.48m), Rectal/Esophageal 24pcs/pk,25°C,10KΩ                                                |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0004FS    | 9FR       |                                                                                                   |      |
|                                                                                                                                                                        |                                                                                     |                                                                                     |         | W0004LS    | 12FR      |                                                                                                   |      |

| OEM Manufacturer & Models                                                                                                                                      | Instrument Connector                                                                | Probe Type                                                                          | OEM #         | Order Code | Dimension | Descriptions                                                                                      | Pice |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------|------------|-----------|---------------------------------------------------------------------------------------------------|------|
| Compatible Mindray (China)                                                                                                                                     |                                                                                     |                                                                                     |               |            |           |                                                                                                   |      |
| PM6800 Series                                                                                                                                                  |    |    | 0011-30-37391 | W0013ES    | 28.8*30mm | Skin-surface Probes<br>24pcs/pk, 25°C, 2.252KΩ                                                    |      |
| BENEVISION/T1/T5/T8<br>IMEC series<br>IPM Series/9800                                                                                                          |                                                                                     |    |               | W0013FS    | 9FR       | 19in. (0.48m),<br>Esophageal/Rectal Probes<br>24pcs/pk, Adult/<br>Pediatric,<br>25°C, 2.252KΩ     |      |
| Compatible Philips                                                                                                                                             |                                                                                     |                                                                                     |               |            |           |                                                                                                   |      |
| All Philips temperature monitors that accept two-prong plugs                                                                                                   |    |    | 21091A        | W0003ES    | 28.8*30mm | 32in.(0.8m),Skin-surface Probes 24pcs/pk,25°C, 2.252KΩ                                            |      |
|                                                                                                                                                                |                                                                                     |    | 21090A        | W0003PS    | 7FR       | 19in.(0.48m),Rectal/<br>Esophageal 24pcs/pk,25°C, 2.252KΩ                                         |      |
|                                                                                                                                                                |                                                                                     |                                                                                     |               | W0003FS    | 9FR       |                                                                                                   |      |
|                                                                                                                                                                |                                                                                     |                                                                                     |               | W0003LS    | 12FR      |                                                                                                   |      |
| Compatible COMEN                                                                                                                                               |                                                                                     |                                                                                     |               |            |           |                                                                                                   |      |
| NC series<br>(NC3/5/8/10/12/19)<br>STARseries<br>(STAR5000A/B/C/D/E/F/H )<br>STAR8000<br>(STAR8000A/B/C/D/E/F/H)<br>C series<br>(C20/30/50/60/70/80/86/90/100) |  |  |               | W0108FS    | 9FR       | 24in. (0.6m),<br>Esophageal/RectalProbes<br>Adult/Pediatric,<br>24pcs/pk,<br>25°C, 10KΩ           |      |
| Compatible Med-Linklet                                                                                                                                         |                                                                                     |                                                                                     |               |            |           |                                                                                                   |      |
| AM-807-C<br>AM-807B<br>AM-807B-C                                                                                                                               |  |  |               | W0024ES    | 28.8*30mm | 32in. (0.8m),<br>Skin-surface Probes<br>24pcs/pk,Adult/Pediatric/Infant/Neonatal<br>25°C, 2.252KΩ |      |
|                                                                                                                                                                |  |  |               | W0114ES    | 28.8*30mm | 32in. (0.8m),<br>Skin-surface Probes<br>24pcs/pk,Adult/Pediatric/Infant/Neonatal<br>25°C, 2.252KΩ |      |
| AM-807<br>AM-808                                                                                                                                               |  |  |               | W0024FS    | 9FR       | 19in. (0.48m),<br>Esophageal/Rectal Probes<br>24pcs/pk, Adult/<br>Pediatric,<br>25°C, 2.252KΩ     |      |

**Medlinket Provide Sterile Disposable Temperature Probes.**