



Water Traps

|
Consumables and
Accessories

Water Traps



CCA003



CCC006



CCB006



CCC007



CCB007

Order Code	Description	OEM #	Package	Scenarios
CCA003	Mini D-Fend	8002174	20/box	ICU/ANES
CCC006	D-Fend	876446-HEL	10/box	ANES
CCB006	D-Fend+	88139-HE	10/box	ICU
CCC007	D-Fend Pro	M1182629	10/box	ANES
CCB007	D-Fend Pro+	M1200227	10/box	ICU

We recommend replacing water traps every 2 months in anesthesia departments and every 24 hours in ICUs. The recommendations are for clinical reference only, and actual replacement intervals should be adjusted based on the infection control risk level and usage conditions of each hospital or department.

Accessories

Product Picture	Order Code	Product Picture	Order Code
	CBE50112		CBB10132
	CBG00112		CBD10112



Shenzhen Med-link Electronics Tech Co., Ltd

Address: Zone A of 1st and 2nd Floor, and 3rd Floor, Building A, No. 7, Tongsheng Industrial Park Road, Shanghenglang Community, Dalang Street, Longhua District, 518109 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
 86-400-058-0755 / 86-0755-61120299
 sales@med-linket.com
 www.med-linket.com

Declaration: Shenzhen Med-Link Electronics Tech Co., Ltd. ("MedLinket") hereby declares that all registered trademarks, product names, model numbers, etc., mentioned in this document belong to respective original owners or manufacturers. The information referenced herein is intended to clearly describe the compatibility features of MedLinket products for reference only and should not be construed as an operational guideline for medical institutions or related entities. If any adverse consequences cause from misunderstandings or incorrect use of the information, MedLinket shall not be held liable.