


Medlinket Products



Holter Recorder ECG Cables with Leadwires

Product Name: Holter Recorder ECG Cables with Leadwires

Brand: Medlinket

Applicable Person: Adult, Pediatric and neonate patients

OEM/ODM: Yes

Quality Certification: CE, FDA, NMPA

*For more product details, check out the information below or contact us directly, marketing@med-linket.com

Compatibility Information:

Compatible Brand	Original Model
Mortara	9293-017-50
GE Healthcare	2008594-001, 2008594-002, 2008594-004, 2008594-002
Philips	M4725A, 350-0173-03, 989803157491
Welch Allyn	SE-PCIEC-PUSH, RE-PCIEC-BAN, RE-PCAHA-BAN, SE-PCIEC-CLIP, SE-PCAHA-CLIP
GE-Marquette Rozinn	E9002TE, E9002TK, E9002TL, E9002TP
Fukuda Denshi	/
Jincomed	/
Creative	/

Specification:


- 1) Leads: 3LD, 4LD, 5LD, 6LD, 7LD, 10LD
- 2) Standard: AHA, IEC
- 3) Patient end terminal: Snap, Clip, Banana



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
 Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
 Email: sales@med-linket.com | Website: www.med-linket.com

Direct-Connect Holter ECG Cables:



Holter ECG Trunk Cables:



ECG Leadwires:



Product Advantages:

1. Trunk Cables fully shielded , Low noise signal transfer;
2. Lighter weight and softer lead wire;
3. Fully compatible with OEM technology;
4. Latex free;
5. Good biocompatibility, free from biological hazard to skin.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of SpO₂, temperature, EEG, ECG, blood pressure, EtCO₂, high-frequency electrosurgical products, etc. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Holter Recorder ECG Cables with Leadwires, ECG Cables with Leadwires, manufacturer, factory

Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com