



4 Infant Incubators

Warmer Temperature Probes

Type Designation

W

Probe type:
C: Reusable Skin-Surface Probes
E: Disposable Skin-surface Probes

Compatible Brand code as 0001
indicated YSI 400 series

Temperature Probes

Probe Type and Code



C

Reusable Skin-surface Probes for Pediatric




E

Disposable Skin-surface Probes


Product Features

1. The accuracy of the thermister is $\pm 0.1^{\circ}\text{C}$ from 25°C to 45°C .
2. Hydrogel stickers available, which is more friendly to babies, and biocompatibility passed.
3. Electromagnetic compatibility.
4. The reflector can effectively insulate ambient temperature and radiant light (Skin-surface).
5. Fully compatible with most of patient monitors.

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Descriptions	Pice
Compatible ATOM (JAPAN)						
Infant Incubator: V-2100, V-2100G, V-2200 Infant Radiant Warmer: V3200, V-3200D/V-505, V-505ST Transport Incubator: V-80, V-800, V-707/V-808, V-850			NF-024/ 60881/ 60881/ 60883	W0138C	118in.(1.5m)Reusable Skin-surface Probes 1pcs/bag	
				W0106C	60in.(3.0m)Reusable Skin-surface Probes 1pcs/bag	
			NF-023/ 60884	W0111E	47in. (1.2m), disposable Skin-surface Probes 24pcs/pk	
Compatible Drager Air-Shields						
ISOLETTE® C2000/8000			MU12525	W0045E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
			MU12525	W0045ES	Sterilized, 47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
			MU12520	W0045K	26in.(0.67m) Adapter cable,5 pin to 4 pin input connection	
			MU12533	W0045C	47in.(1.2m), reusable Skin-surface Probes,1pcs/bag	
IsoletteC450 C-550 QT, Isolette C100, Isolette C100 QT, Isolette C450/QT, Isolette C550, GLOBE-TROT_x005f TER GT5400, IICS-90, PM78			MU06951	W0010C	60in.(3.0m), reusable Skin-surface Probes,1pcs/bag	
			MU06933	W0010E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
Incubators (BABYTHER_x005f M®)8000 WB/8004/8010/ 4200 CALEO®				W0020C	32in.(0.8m), Reusable Skin-surface Probes,1pcs/bag	
			MX11001	W0020E	32in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
			MX11001	W0020ES	Sterilized, 32in. (0.8m), disposable Skin-surface Probes, 24pcs/pk	
				W0095E	63in.(1.6m), disposable Skin-surface Probes,24pcs/pk	

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Descriptions	Pice
Compatible David (China)						
InfantIncubator: YP-3000, YP-2008, YP-2000, YP-970, YP-930/920, 910 (The above is a new machine)			88.10.005.00	W0119C	59in.(1.5m) reusable Skin-surface Probes,1pcs/bag	
				W0119E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
			88.09.005.00	W0120C	59in.(1.5m), reusable Skin-surface Probes,1pcs/bag	
				W0120E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
Infant Radiant Warmer: HKN-93/93A/93B/ 93C/, HKN-2001, HKN-9010			88.10.005.00	W0119C	59in.(1.5m) reusable Skin-surface Probes,1pcs/bag	
				W0119E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
Transport Incubator: TI-2000			88.09.005.00	W0120C	59in.(1.5m), reusable Skin-surface Probes,1pcs/bag	
				W0120E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
Compatible GE Datex-ohmeda						
Giraffe series			6600-0875-700	W0026C	118in.(3.0m), reusable Skin-surface Probes,1pcs/bag	
			6600-0873-700 6600-0874-700	W0026E	118in.(3.0m), disposable Skin-surface Probes,24pcs/pk	
General Care (GC) Incubators (3000/3300/4000/5000), Ohio Infant Care (IC), Panda Infant Warmers			0208-0697-700 6600-0628-700	W0021C	59in.(3.0m) reusable Skin-surface Probes,1pcs/bag	
			6600-208-700 6600-0196-700	W0021E	47in.(3.0m) disposable Skin-surface Probes,24pcs/pk	

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Descriptions	Pice
Compatible Zhengzhou Dison						
Infant Incubator: BB-300W/300G/ 300/200W/200G/ 200/100W/100G/ 100 Infant Radiant Warmer: BN-200/100/100A				W0117C	59in.(1.5m), reusable Skin-surface Probes,1pcs/bag	
				W0117E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
				W0118C	59in.(1.5m), reusable Skin-surface Probes,1pcs/bag	
				W0118E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
Compatible Julongsanyou						
Infant Incubator: B-4500/4000/ 3000/2000/ 1000/800				W0123C	59in.(1.5m), reusable Skin-surface Probes,1pcs/bag	
				W0123E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	
Infant Radiant Warmer: BNT-300/2000B/ 2000/1000				W0122C	59in.(1.5m), reusable Skin-surface Probes,1pcs/bag	
				W0122E	47in.(1.2m), disposable Skin-surface Probes,24pcs/pk	

Temperature Probe Accessories				
Order Code	Product Picture	Descriptions	Scope of Application/ Characteristics	Pice
BZ0109D		Heart shape 28.8*30mm, reflective radiating patch	The radiation reflection paste is used in combination with reusable and disposable temperature probe, which is fixed on human epidermis to ensure more accurate body temperature detection. Characteristic: 1) Which can effectively isolate ambient temperature and radiated light; 2) Optional hydrogel radiation reflective patches that are more friendly to infants.	
BZ0501A		Circular shape, 30mm, reflective radiating patch		
BZ0109F		Heart shape, 28.8*30mm, hydrogel reflective radiating patch		