

## Air Hose Connectors (Cuff Side)



**Product Name:** Air Hose Connectors

**Brand:** Medlinket

**OEM/ODM:** Yes

**Descriptions:** Fit for ID 4.0 mm tube, Fit for ID 2.5-4.0 mm tube

\*For more product details, check out the information below or contact us directly

### Air Hose Connectors (Cuff Side):

1. Air Hose Connectors (Suitable for Adult / Pediatric):

Compatible Brand	OEM#
Air Hose Connector	GE P/N: 300670
L&T, Critikon, Datex	GE P/N: 300669 Welch-allyn P/N; 5082-169
GE, Marquette and Protocol	GE P/N: 330064 Welch-allyn P/N; 5082-182
Philips, Siemens, Datascope and Colin	GE P/N: 330059, 330060; Welch-allyn P/N: 5082-184
Spacelabs, Datascope and Colin	GE P/N: 300668, 300665; Welch-allyn P/N; 5082-168, 5082-165
GE DINAMAP	GE P/N: 300664, 300619; Welch-allyn P/N; 5082-164, 5082-161
Nihon Kohden SVM/BSM/PVM Series	/
GE Marquette Datex-Ohmeda	GE P/N: 330090
GE Marquette one-tube systems	GE P/N: 330057; Welch-allyn P/N: 5082-181
Welch Allyn	16-41-000, 16-40-000
ADAPTER	/
GE	212202296
Drager	/



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
Email: sales@med-linket.com | Website: www.med-linket.com



## 2. Air Hose Connectors (Suitable for Infants / Newborns)

- 1) Compatible Brand: GE, Datex, Welch-allyn, Pacific Medical, Philips
- 2) OEM #: 300667, 5082-176, PM20



### Product Advantages:

1. For adults, children, infants and newborns;
2. Material hose to ensure good air tightness and prolong service life;
3. Fully compatible with original monitor blood pressure connector;
4. Good biocompatibility, free from biological hazard to skin;
5. Latex free, PVC free.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of air hose connectors in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: [marketing@med-linket.com](mailto:marketing@med-linket.com).

**Hot Tags:** Air Hose Connectors, Air Tubing, Y Type Tube Adapter

\*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
Email: [sales@med-linket.com](mailto:sales@med-linket.com) | Website: [www.med-linket.com](http://www.med-linket.com)