


Disposable Anesthesia EEG Sensor



Product Name: Disposable Anesthesia EEG Sensor

Brand: Medlinket

Applicable Person: Adult, Pediatric

Sensor Type: 3 electrode, 4 electrode, 6 electrode

OEM/ODM: Yes

Quality Certification: CE, NMPA, ISO9001, ISO13485

*For more product details, check out the information below or contact us directly.

Product Advantages :

1. No exfoliation with sandpaper, small workload and avoiding resistance detection failure due to improper wiping;
2. Small electrode size, no effect for the attachment of the brain oxygen probe;



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

3. Single patient use, avoiding cross-infection;
4. High-quality conductive plastics and sensor parts, quick response;
5. Good biocompatibility, free from biological hazard to skin;
6. Optional waterproof sticker device.

Compatibility Information:

OEM #	Adaptable Module
186-0106	Double channel module Bispectral index (Adult)
186-0200	Double channel module Bispectral index (Pediatric)
/	EEG state index
M1174413	Entropy index
2479	Compatible Masimo SedLine® Brain Function Monitoring For use with SedLine MOC-9 and and SedLine Patient Cable
186-0212	Four channel module Bispectral index
/	IOC anesthesia depth

Anesthesia Adapters:



As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of disposable anesthesia EEG sensor in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Disposable Anesthesia EEG Sensor

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com