

## Disposable NIBP Cuff



**Product Name:** Disposable NIBP Cuff

**Brand:** Medlinket

**OEM/ODM:** Yes

\*For more product details, check out the information below or contact us directly

### Product Advantages:

1. For single patient only, avoiding cross infection;
2. Soft and comfortable, Low hazard to skin even if a long-term use;
3. Transparent design for observing the patient skin condition conveniently;
4. Easy to use, The index line for proper size and placement;
5. Latex free, DEHP free;
6. Good biocompatibility, free from biological hazard to skin.

### Scope of Application:

1. Operating Room (OR), ICU, Neonatology Internal Cardiovascular Dept, Cardiothoracic Surgery Dept;
2. For patients prone to infection such as burns, open surgery, newborns, infectious diseases, etc.

### Disposable Neonate NIBP Cuff:

(1) Disposable NIBP Soft Fiber Cuff/Hylink Disposable NIBP Comfort Cuff

Limb Circumference	Single Tube	Double Tube
	OEM #	OEM #
3-6 cm	5082-101-1	5082-101-2
4-8 cm	5082-102-1	5082-102-2
6-11 cm	5082-103-1	5082-103-2
7-14 cm	5082-104-1	5082-104-2



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
Email: sales@med-linket.com | Website: www.med-linket.com

8-15 cm	5082-105-1	5082-105-2
---------	------------	------------

## (2) Compatible Philips Disposable NIBP Comfort Cuff

Limb Circumference	Single Tube
	OEM #
3-6 cm	M1866B
4-8 cm	M1868B
6-11 cm	M1870B
7-14 cm	M1872B
8-15 cm	M1873B

## Disposable NIBP Cuff:

Disposable NIBP Comfort Cuff:

Patient Size	Limb Circumference	Single Tube	Double Tube
		OEM #	OEM #
Adult thigh	42-50 cm	5082-98-3	5082-98-4
Large adult	32-42 cm	5082-97-3	5082-97-4
Adult long	28-37 cm	5082-96L-3	5082-96L-4
Adult	24-32 cm	5082-96-3	5082-96-4
Small adult	17-25 cm	5082-95-3	5082-95-4
Pediatric	15-22 cm	5082-94-3	5082-94-4

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of disposable NIBP cuff in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: [marketing@med-linket.com](mailto:marketing@med-linket.com).

Hot Tags: Disposable NIBP Cuff

\*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
Email: [sales@med-linket.com](mailto:sales@med-linket.com) | Website: [www.med-linket.com](http://www.med-linket.com)