


ECG Cable and Leadwire (for OR)



Product Name: ECG Cable and Leadwire (for OR)

Brand: Medlinket

OEM/ODM: Yes

Compatible Brand: Philips

OEM #: 989803170171, 989803170181

*For more product details, check out the information below or contact us directly

Product Advantage:

★Good anti-electric knife performance, used in operating room, electric cut signal quickly recovery;


★The cable connect with it can replaced,convenient and flexible;

★Used soft and comfortable TPU cable,can be repeatedly cleaned and disinfected;

★Good biocompatibility, latex free, avoid allergic reactions to patients.

Description:

1) ECG trunk cable; 2) 12pin; 3) 3-Lead ESU-Proof, 5-Lead ESU-Proof.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

Scope of Application:

Connected between instrument and electrode, used with ECG monitor and electrocardiogram electrode, to transmit electrophysiological signals collected from the body surface.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of ECG cable and leadwire (for OR) in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: ECG Cable and Leadwire (for OR), ECG Trunk Cable

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com