


ESU Pencil



Product Name: ESU Pencil

Brand: Medlinket

OEM/ODM: Yes

Compatible Brand: Erbe, Valleylab, Bovie, BOWA, Conmed, Martin

*For more product details, check out the information below or contact us directly

Descriptions:

- 1) Disposable Electrosurgical Pencil, Reusable Electrosurgical Pencil
- 2) Blade Length=40mm
- 3) Length of wire: 2.8m, 3m, 5m
- 4) Cable material: PVC, silicone
- 5) Plug: 3pin banana plug, $\phi 4.83$ banana plug (gold-plated), $\phi 4.0$ banana plug
- 6) sterilized by ethylene oxide

Product Advantages:

1. Ergonomic design, comfortable feel and excellent mechanical properties;
2. Different shapes of cutting heads can be changed according to the needs of surgery;
3. Closed structure and unique anti-rotation hexagon design, safe and reliable;
4. Extended and increased anti-folding design, strong flexural folding ability, easy to clean, multiple protection;
5. Adaptable plugs can be customized according to adapted model;
6. V type design, easy to plug and unplug, reduced the risk of surgical maximum;
7. Silicone wire, steam sterilizable.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com

We can supply our customer with common active electrode and customize special active electrode according to their requirements.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of ESU pencil in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: ESU Pencil, Disposable Electrosurgical Pencil, Reusable Electrosurgical Pencil, ESU

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com