


Electrosurgical Device Cable



Product Name: Electrosurgical Device Cable

Brand: Medlinket

OEM/ODM: Yes

Compatible Brand: Erbe, T-seris, Martin, Berchtold, Wolf, Bovie, Valleylab, Conmed, ACC/ICC, Gyrus Acmi

OEM#: 20192-117, 815.033, 8106.034, 3900

*For more product details, check out the information below or contact us directly

Descriptions:

- 1) Length of wire: 9.8ft (3m), 10ft (3m), 13ft (4m), 15ft (4.5m), 16ft (5m)
- 2) Cable material: silicone
- 3) Connector type: 7.9 mono plug, 8.0 mono plug, $\Phi 4.83$ banana plug, $\Phi 4.0$ copper tube female plug, $\Phi 4.0$ assembled banana plug, $\Phi 4.0$ red head banana female plug.

Instrument Connector:



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

Patient Connector:



As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of electrosurgical device cable in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Electrosurgical Device Cable, ESU

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com