

## EtCO2 Mainstream Module



**Product Name:** EtCO2 Mainstream Module

**Brand:** Medlinket

**Applicable Person:** Adult, Pediatric, Infant

**Compatible Brand:** Respironics, Masimo, ZOLL (E/R Series), Philips, Mindray (China)

**OEM #:** 1015928, 200601(IRMA AX+), 8000-0312, M2501A 989803142651, 6800-30-50760

**Descriptions:** 8pin, 9pin; Mainstream EtCO2 Sensor

\*For more product details, check out the information below or contact us directly

### Product Introduction:

Medlinket provides a cost-effective EtCO2 monitoring program for clinical practice. It is plug-in and play. Advanced non-spectrophotometric infrared technology can be used to measure the instantaneous CO2 concentration, respiratory rate, end-tidal CO2 value and inspired CO2 concentration of the measured object.

### Connector Type:



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

**Product Features:**

1. Simple Operation;
2. Stable, dual a1 waveband, NDIR (non-dispersive infrared) technology;
3. Long service life, infrared biakbody light source of MEMS technology;
4. Accurate calculation result, compensating temperature, pressure and Bayesian gas;
5. Calibration-free, patented calibration algorithms;
6. With a minimum samplingflow rate of 50ml/min;
7. Srong compatibility, adapt for various brand modules.

**Application Field:**

1. Monitoring the patient's respiratory status;
2. Assists in determining when to intubate or extubate;
3. Verification of ET tube placement;
4. Alerts if accidental extubation occurs;
5. Ventilator disconnect detection;
6. Verification of ventilation during transport.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of reusable spo2 sensor in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: [marketing@med-linket.com](mailto:marketing@med-linket.com).

Hot Tags: EtCO2 Mainstream Capnography, EtCO2 Sidestream Module, EtCO2 Mainstream Capnography

\*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



**Shenzhen Med-link Electronics Tech Co., Ltd**

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055  
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China  
Email: [sales@med-linket.com](mailto:sales@med-linket.com) | Website: [www.med-linket.com](http://www.med-linket.com)