

## Flow Sensor Cable



**Product Name:** Flow Sensor Cable

**Brand:** Medlinket

**OEM/ODM:** Yes

**Compatible brand:** Draeger Babylog Vn500 ventilator

**OEM#:** 8416600

**Specification:** 7Pin>4Pin , length 4ft

\*For more product details, check out the information below or contact us directly

### Product Advantage:

- ★Net tail anti-dusty design,easy to clean;
- ★Plug connector ergonomic design,easy to plug and unplug;
- ★Good biocompatibility, latex free, avoid allergic reactions to patients;

### Scope of Application:

Used to connect flow sensor (P/N:8403735,6871980,MK01900,8410185,8411130,8410179)

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of flow sensor cable in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: [marketing@med-linket.com](mailto:marketing@med-linket.com).

Hot Tags: Flow Sensor Cable, Flow sensor

\*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



### Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: [sales@med-linket.com](mailto:sales@med-linket.com) | Website: [www.med-linket.com](http://www.med-linket.com)