

Inspiratory Heater Wire Adapter



Product Name: Inspiratory Heater Wire Adapter

Brand: Medlinket

OEM/ODM: Yes

*For more product details, check out the information below or contact us directly

Product Advantage:

- ★ Ergonomically design, easier plug-in connection.;
- ★ Latex free, cost-effective.

Compatibility Information:

Compatible Brand	Fisher & Paykel 700 series, HC550 humidifier, MR850 humidifier, RT-Series Breathing Circuit, House heater 900MR755(20Ω)9.4Ω/m
#OEM	900MR858
Specification	Length 1.6ft (0.5M) 4-pin connector >> 2-pin connector
Related Products	W0133Y, W0133M, W0127I

Scope of Application:

Connected with respiratory humidifier and heated breathing circuit and used for heating of single breathing circuit to avoid condensate water in the tube line.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of inspiratory heater wire adapter in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com

certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Inspiratory Heater Wire Adapter

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com