


NIBP Adapter Air Hose



Product Name: NIBP Adapter Air Hose

Brand: Medlinket

OEM/ODM: Yes

Applicable Person: Adult, Pediatric, Infant, Neonate

*For more product details, check out the information below or contact us directly

Descriptions:

1) Length: 1.5m, 1.8m, 2m, 2.4m, 3m, 3.6m, 7.3m

2) Diameter: 4 mm, 2.5 mm-4.0 mm

Compatibility Information:

Compatible Brand	OEM#
Bionet (Korea)	/
Criticare/CSI	/
Colin	/
Creative (China)	/
Drager/Siemens	MP00953, 2870298, MP00953
Fukuda	/
Goldway (China)	/
GE	107365, 107366, 107368, 107363, 9461-203, 9461-214, 9461-205, 9461-217, 414873-001, 2058203-002, 2017009-001, 414874-001, 88845, 2087389-002
Mindray (China)	6200-30-09688, 6200-30-11560,
Nihon Kohden	YN-900P, YN-920P, YN-701S
Neusoft Jeteem	/
Medtronic Physio Control	/
Philips	M1598B, M1599B, M1597C, M3918A



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

Unicare Medical Centre	/
Welch Allyn	4500-30, 4500-31
Siemens	/
ZOLL	8300-0002-01, 8000-0655
Spacelabs	714-0018-00

Plug to cuff:



Plug to monitor:



Product Advantages:

1. It will not be deformed even if long-term use;
2. Reasonable hardness avoids severe bending, ensures the safety of deflation;
3. Reasonable hardness ensures well transmission of blood pressure signal;
4. TPU material hose to ensures good air tightness and long life;
5. Latex free, PVC free;
6. Good biocompatibility, free from biological hazard to skin.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of NIBP adapter air hose in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: NIBP Adapter Air Hose

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com