

Pressure Infusion Bags



Product Name: Reusable Pressure Infusion Bags, Disposable Pressure Infusion Bags

Brand: Medlinket


OEM/ODM: Yes

Material: nylon, waterproof

OEM #: MX4805, MX4810, MX4830

*For more product details, check out the information below or contact us directly.

Product Features:



Descriptions:

- 1) Capacity: 500ml, 1000ml, 3000ml
- 2) L*W: 309*150mm, 380*150mm, 380*220mm
- 3) Silicone tube, PVC pipe



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

Scope of Application:

It is suitable for emergency blood transfusion and fluid infusion in various clinical departments such as emergency dept., operating room, ICU, emergency center, etc.

Product Advantages:

1. Reuse by single patient to reduce cross-infection (disposable);
2. Medical PU composite material, comfortable to use and easier to clean (disposable);
3. Waterproof nylon cloth, not easy to tear; silicone inflatable tube, comfortable and easy to use (reusable);
4. Pressure indicator, precise control of pressure and flow;
5. Balloons with soft texture, good flexible, palm-sized, and efficient inflation;
6. Over-pressure protection device to avoid excessive inflation pressure and burst, frightening the patient;
7. Equipped with Robert clamp, more reliable in pressure holding;
8. The unique hook design, avoiding the risk of falling off after the blood bag or liquid bag volume is reduced, and assuring safer to use.

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of pressure infusion bags in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Pressure Infusion Bags, Reusable Pressure Infusion Bags, Disposable Pressure Infusion Bags

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com