

Reusable Temperature Probes



Product Name: Reusable Temperature Probes

Brand: Medlinket

Applicable Person: Adult, Pediatric, Neonatal

Resister Type: 25 °C 6.867 K Ω , 25 °C 2.252 K Ω , 25 °C 10 K Ω , 25 °C 2.016K Ω , 18B20

Descriptions: 10 ft (3.0 m)

*For more product details, check out the information below or contact us directly

Compatible Brand: Air Shields; Artema/S&W, Comen (China); Creative (China); Drager; Edan (China); Mindray-Datascope; Mindray (China); Philips; Spacelabs; YSI 400 Series; 10 k, 25 °C; 2.016 K, 25°C; DALLAS 18B20 Probe

OEM #: T1302, 040-000246-00, 040-000385-00, 040-000387-00, 040-000386-00, 5204669, 4329889, 01.15.040225, 01.15.040228, 01.15.040255, 01.15.0402256, 040-000057-00, 040-000055-00, 040-000058-00, 040-000056-00, 0011-30-37393, (0011-30-37382)(0011-30-37392), 0011-30-37395(MR404B), 0011-30-37394, 21078A, 21075A, 21076A, 20700-4000-00A, 409B,709B, 701,85409, YSI: 401 GE P/N: 9445-004, YSI: 427 Welch Allyn P/N: 008-0756-00, YSI: 402 GE P/N: M1024251, Ref: 092.059.600

Probe Type: Reusable Skin-surface Probes, Reusable Esophageal/Rectal Probe



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

Instrument Connector:



Product Advantages:

1. High precision: ± 0.1 °C
2. Quick response
3. High performance noise ratio
4. Fully compatible with OEM technology
5. Complete product specifications

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of Reusable Temperature Probes in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Medical Temp Probe, Monitor temperature probe, Reusable Temperature Probe, YSI 400 series Temp

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above. information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055
Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China
Email: sales@med-linket.com | Website: www.med-linket.com