


Sterile Disposable Subdermal Needle Electrodes



Product Name: Sterile Disposable Subdermal Needle Electrodes

Brand: Medlinket

OEM/ODM: Yes

*For more product details, check out the information below or contact us directly.

Product Advantages:

1. Sterilized by EO gas, Single patient use only;
2. To provide easy and secure recording when placed directly into the scalp;
3. Made of stainless steel and designed for strength and flexibility;
4. The tip is sharpened with precision to an ultra-sharp lancet cut for easy penetration;
5. The twisted lead wires reduce noise and impedance levels giving a clean, clear signal.

Descriptions:

1. Electrode length: 12mm/0.5in, 14mm/0.55in
2. Electrode diameter: 0.40mm/27gauge
3. Plug, Leadwire:
 - 1) Single wire, Twisted pairs
 - 2) DIN Style, DIN 1.5
 - 3) 1.5m
 - 4) Color: green, red, white, purple, blue, yellow, brown, orange, black



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com

As a professional manufacturer of various quality medical sensors & cable assemblies, Med-linket is also one of the leading suppliers of sterile disposable subdermal needle electrodes in China. Our factory is equipped with advanced equipment and many professionals. With FDA and CE certification, you can rest assured to buy our products made in China at reasonable price. Also, OEM / ODM customized service are also available.

If you need more information, please feel free to contact us: marketing@med-linket.com.

Hot Tags: Sterile Disposable Subdermal Needle Electrodes, EMG Electrode

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company.



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com