



3 Reusable Temperature Probes

Type Designation

W

Probe Type:

A: Adult reusable skin-surface probes
 B: Adult reusable esophageal/rectal probe
 C: Pediatric reusable skin-surface probes
 D: Pediatric reusable esophageal/rectal probe
 N: Reusable skin surface probes for Neonate

Compatible brand code as 0001
 indicated YSI 400 series

Temperature Probes

Probe Type and Code



Reusable skin surface probes for adult



Reusable esophageal/rectal probe



Reusable skin surface probes for Pediatric



Reusable esophageal/rectal probe for Pediatric



Reusable skin surface probes for Neonate

Product Features

1. High precision: $\pm 0.1^{\circ}\text{C}$
2. Quick response
3. High performance noise ratio
4. Fully compatible with OEM technology
5. Complete product specifications

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Resister Type	Descriptions	Price
Temperature accessories for monitor							
Compatible Air Shields							
C-550 QT, Isolette C100/2000/C450, Isolette C100 QT/C450 QT/C550, PM78				W0010A	25°C 6.867KΩ	10ft (3.0m), Adult skin	B0
				W0010B	25°C 6.867KΩ	10ft (3.0m), Adult Rectal/Esophageal	B0
Compatible Artema/S&W							
DIASCOPE G2			T1302	W0007B	25°C 2.252KΩ	10ft (3.0m), Adult Rectal/Esophageal	B5
Compatible Comen (China)							
STAR 8000 series, NC series, K series, C series				W0092A	25°C 2.252KΩ	10ft (3.0m), Adult skin	B3
				W0092B	25°C 2.252KΩ	10ft (3.0m), Adult Rectal/Esophageal	B3
				W0092C	25°C 2.252KΩ	10ft (3.0m), Pediatric skin	B3
				W0092D	25°C 2.252KΩ	10ft (3.0m), Pediatric Rectal/Esophageal	B3
			040-000246-00	W0128A	25°C 10KΩ	10ft (3.0m), Adult skin	B3
			040-000385-00	W0128B	25°C 10KΩ	10ft (3.0m), Adult Rectal/Esophageal	B3
			040-000387-00	W0128C	25°C 10KΩ	10ft (3.0m), Pediatric skin	B3
			040-000386-00	W0128D	25°C 10KΩ	10ft (3.0m), Pediatric/ Neonatal Rectal/Esophageal	B3
Compatible Creative (China)							
PC-1000, 5000 9000, UP-6000				W0005A	25°C 10KΩ	10ft (3.0m), Adult skin	B5
Compatible Biolight							
Biolight A series (A2/A3/A5/A6/A8), Q series (Q3/Q4/Q5/Q6/Q7)			15-031-0005	W0113A	25°C 2.252KΩ	10ft (3.0m), Adult skin	B3
			15-031-0012	W0113B	25°C 2.252KΩ	10ft (3.0m), Adult Rectal/Esophageal	B3
				W0113C	25°C 2.252KΩ	10ft (3.0m), Pediatric skin	B3
				W0113D	25°C 2.252KΩ	10ft (3.0m), Pediatric Rectal/Esophageal	B3

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Resister Type	Descriptions	Price
Temperature accessories for monitor							
Compatible Drager							
Infinity Delta, Infinity Gamma, Infinity Kappa, Infinity Vista			5204669	W0008A	25℃ 2.252KQ	10ft (3.0m),Adult skin	B5
			4329889	W0008B	25℃ 2.252KQ	10ft (3.0m),Adult Rectal/Esophageal	B5
				W0008C	25℃ 2.252KQ	10ft (3.0m), Pediatric skin	B7
				W0008D	25℃ 2.252KQ	10ft (3.0m),Pediatric Rectal/Esophageal	B7
Compatible Edan (China)							
IM series (50/60/70/80); Elite V series (V5/V5/V8); IX series(X8/X10/ X12/X15)			01.15.0402 25	W0116A	25℃ 10KQ	10ft (3.0m),Adult skin	C0
			01.15.0402 28	W0116B	25℃ 10KQ	10ft (3.0m),Adult Rectal/Esophageal	C0
			01.15.0402 55	W0116C	25℃ 10KQ	10ft (3.0m), Pediatric skin	C0
			01.15.0402 256	W0116D	25℃ 10KQ	10 ft(3.0m),Pediatric Rectal/Esophageal	C0
Compatible GE							
GE CARESCAPE series (B450/B650/ B850), DASH series (2500/3000/4000), B20/40series, B105M/B105P/ B125M/B125P/ B155M				W0015A	25℃ 2.252KQ	10ft (3.0m),Adult skin	B0
				W0015B	25℃ 2.252KQ	10ft (3.0m),Adult Rectal/Esophageal	B0
				W0015C	25℃ 2.252KQ	10ft (3.0m), Pediatric skin	C1
				W0015D	25℃ 2.252KQ	10ft (3.0m),Pediatric Rectal/Esophageal	C1
Compatible Mindray-datascope							
Passport V, V-Series			040-0000 57-00	W0112A	25℃ 2.252KQ	10ft (3.0m),Adult skin	C0
			040-0000 55-00	W0112B	25℃ 2.252KQ	10ft (3.0m),Adult Rectal/Esophageal	C0
			040-0000 58-00	W0112C	25℃ 2.252KQ	10ft (3.0m), Pediatric skin	C0
			040-0000 56-00	W0112D	25℃ 2.252KQ	10ft (3.0m),Pediatric Rectal/Esophageal	C0

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Resister Type	Descriptions	Price
Temperature accessories for monitor							
Compatible Mindray (China)							
BeneVision N Series: N1/N17/N15/N12/N2 2/N19/N15/N17 ICU; BeneVision(BeneV-iew) T Series: T1/T5/T6/T8/T9; ePM Series: ePM10/12/15/10M/ 12M;iPM Series: iPM8/10/12; iMEC Series:iMEC8/10/12/5C ;uMEC Series: uMEC10/12			0011-30-37393	W0013A	25℃ 2.252KQ	10ft (3.0m),Adult skin	B0
			(0011-30-37382) (0011-30-37392)	W0013B	25℃ 2.252KQ	10ft (3.0m),Adult Rectal/Esophageal	B0
			0011-30-37395 (MR404B)	W0013C	25℃ 2.252KQ	10ft (3.0m), Pediatric skin	C1
			0011-30-37394	W0013D	25℃ 2.252KQ	10ft (3.0m),Pediatric Rectal/Esophageal	C1
Compatible Philips							
Philips: MP Series: MP20/30/40/50/60/70 /80/90; MX Series: MX400/450/500/550; SureSigns VM4/VM6/VM8; Goldway: CM Series: CM100/10/120/12/ 150			21078A	W0003A	25℃ 2.252KQ	10ft (3.0m),Adult skin	B5
			21075A	W0003B	25℃ 2.252KQ	10ft (3.0m),Adult Rectal/Esophageal	B5
				W0003C	25℃ 2.252KQ	10ft (3.0m), Pediatric skin	B7
			21076A	W0003D	25℃ 2.252KQ	10ft (3.0m),Pediatric Rectal/Esophageal	B7
Compatible Sunray							
Sunray: SRF				W0134A	25℃ 2.252KQ	10ft (3.0m), Adult skin	B0
Compatible Spacelabs							
Spacelab: New modules 90651-08, 90496, Ultraview			20700-4000-00A	W0022A	25℃ 2.252KQ	10ft (3.0m), Adult skin	B0
				W0022B	25℃ 2.252KQ	10ft (3.0m), Adult Rectal/Esophageal	B0
Compatible YSI 400 series							
CAS Med: 750 CSI: 1100, 2200 Ohmeda: CS/3, Cardiicap I Draeger: Infinity Delta, Infinity Gamma Edan: M50, M80 Datascope: PM 7000 Nihon Kohden: Spacelabs: Welch Allyn: YSI: 400 Series Zoll			409B, 709B, 701, 85409	W0001A	25℃ 2.252KQ	10ft (3.0m), Adult skin	A9
			YSI: 401 GE P/N: 9445-004	W0001B	25℃ 2.252KQ	10ft (3.0m), Adult Rectal/Esophageal	A9

OEM Manufacturer & Models	Instrument Connector	Probe Type	OEM #	Order Code	Resister Type	Descriptions	Price
Temperature accessories for monitor							
Compatible YSI 400 series							
CAS Med: 750, CSI: 1100, 2200, Ohmeda: CS/3, Cardiocap I, Draeger: Infinity Delta, Infinity Gamma, Edan: M50, M80 Datascope: PM 7000 Nihon Kohden: Spacelabs:Welch Allyn: YSI: 400 Series, Zoll			YSt:427 Welch Allyn P/N:008-07 56-00	W0001C	25°C 2.252KQ	10ft (3.0m), Pediatric skin	A0
			YSI:402 GE P/N: M1024251	W0001D	25°C 2.252KQ	10ft(3.0m),Pediatric Rectal/Esophageal	A0
				W0029A	25°C 2.252KQ	10ft(3.0m),Adult skin	A9
				W0029C	25°C 2.252KQ	10ft (3.0m), Pediatric skin	A0
YSI 400 series				W0029N	25°C 2.252KQ	10ft (3.0m),Neonate 6.35mm Mono Plug	A0
YSI 400 (Need Extension)				W0001N	25°C 2.252KQ	5ft (1.5m), Neonate Skin-surface Probes	B7
Compatible 10k, 25°C							
Tianrong: TR900 D/E; Anke: ASC553A3, ASC553; Comen: star 8000A/B/C, star 5000, 5000B/C				W0004A	25°C 10KQ	10ft (3.0m),Adult skin	A8
				W0004B	25°C 10KQ	10ft (3.0m),Adult Rectal/Esophageal	A8
Compatible 2.016K, 25°C							
YSI 25°C@2.016KQ O.S.7929			Ref.092.0 59.600	W0009A	25°C 2.016KQ	10ft (3.0m),Adult skin	B0
Compatible DALLAS 18B20 Probe							
18B20 series				W0006A	18B20	10ft (3.0m),Adult skin	B0
				W0006B	18B20	10ft (3.0m),Adult Rectal/ Esophageal	B0
Compatible Med-linket							
AM-801, AM-805, AM-806, AM-806B				W0114C	25°C 2.252KQ	3.3ft(1.0m),Pediatric Rectal/Esophageal	B7
				W0001N (Need Adapters Cable S0177MD- S)	25°C 2.252KQ	5ft(1.5m), Neonate Skin-surface Probes	B7