

New Product



IBP Adapter Cable(For BD Transducer)

Product Advantage

- ★One-piece injection molding process, strong and durable, easy to clean;
- ★The connector feels good when plugging and unplugging, soft and comfortable;
- ★No latex, cost-effective;
- ★Excellent shielding performance and anti-interference performance, better signal quality.

Scope of Application

After connecting to the multi-parameter monitor, measure the patient's arterial and venous blood pressure.

Product Parameter

Compatible Models	Nihon Kohden Life Scope series; BSM series		
Brand	Medlinket	Order code	X0045D
Specification	5Pin>4Pin,length 12ft(3.6m)	Color	Cold gray
Weight	164g/pcs	Related product	X0046D
Package	1pcs/bag		

*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company



Shenzhen Med-link Electronics Tech Co., Ltd

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor,Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com