

# New Product



## Blood Pressure Cuff

### Product Advantage

- ★The jacket is made of nylon fabric, which is soft, skin-friendly and comfortable;;
- ★TPU inner tank, good air tightness, accurate test and long service life;
- ★Clear range markings and operating steps for easy selection of the appropriate cuff and ease of use;
- ★Compatible with Omron Series 5 electronic sphygmomanometers, a cost-effective alternative;
- ★Good biocompatibility, latex free, avoid allergic reactions to patients.

### Scope of Application

Through vasoconstriction and expansion, the pressure of the cuff liner is collected and transmitted the human blood pressure signal, which is suitable for general wards of hospitals, clinics and households.

### Product Parameter

Compatible Brand	Omron Series 5		
Picture	Order Code	Limb Circumference	Descriptions
A	Y003A1-A62	22-32CM	Suitable for adults, single tube, trachea length: 61.5 cm, nylon fabric.
B	Y003L1-A62	32-45CM	Suitable for adult plus size, single tube, trachea length: 61.5 cm, nylon fabric.

\*Declaration: All the registered trademarks, names, models, etc. displayed in the above content are owned by the original owner or original manufacturer. This article is only used to illustrate the compatibility of the Med-Linket products. There is no other intention! All the above information is for reference only, and should not be used as a guide for the work of medical institutions or related units. Otherwise, any consequences caused by this company have nothing to do with this company



**Shenzhen Med-link Electronics Tech Co., Ltd**

Tel: (86) 755-61120299 | Oversea Dept: (86) 755-61120085 | Fax: (86) 755-61120055

Address : 4th and 5th Floor, Building Two, Hualian Industrial Zone, Xinshi Community, Dalang Street, Longhua District, 518109 Shenzhen, China

Email: sales@med-linket.com | Website: www.med-linket.com